Technical Review

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Feature of Extreme Liquid Chromatography System X-LC™
and Application to Ultra High-Speed Separation of Amino Acids

Masao Bounoshita, Yasuyo Sato, Toshihiko Miyaji, Keijin Iwaya,
Higashidate Sakae, Takayuki Yamaguchi and Muneo Saito

JASCO Corporation, 2967-5 Ishikawa-cho, Hachioji-shi, Tokyo, 192-8537 Japan

Abstract

Ultra high-speed liquid chromatography has become increasingly popular in HPLC. In general, conventional HPLC separation has been performed on columns packed with 5 μm particles. However, conventional HPLC requires a long analysis time and a large amounts of solvents. In case of ultra high-speed liquid chromatography, shorter columns typically 50 – 100 mm in length and 2 mm I.D. packed with ultra micro-particles, of which particle size is less than 2 μm, are used. X-LC (extreme liquid chromatograph) offers higher efficiency and sensitivity, and reduction of organic solvents without sacrificing separation efficiency. We outline the new generation HPLC equipment X-LC that can be used for ultra high-speed liquid chromatographic separation. We also discuss the analysis of amino acids using the X-LC system, utilizing a combination of precolumn derivatization with OPA (orthophthalic aldehyde) and separation on a 2-μm dia C18 column with fluorescence detection.

Keywords: Ultra high-speed liquid chromatography, X-LC, Extreme liquid chromatography, Amino acids, precolumn derivatization, OPA, preheat-coil, High-Throughput

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Eluent A: 1.0 M citrate buffer (pH 5.8) 3.5 mL in 1 L of H₂O

Eluent B: 1.0 M citrate buffer (pH 5.8) 3.5 mL in 1 L of CH₃CN/C₂H₅OH/H₂O (30/30/40)

Gradient condition : 1 cycle 10 min

A:B = 90:10 $\xrightarrow{0.2 \text{ min}}$ 90:10 $\xrightarrow{2.2 \text{ min}}$ 72:28 $\xrightarrow{2.5 \text{ min}}$ 72:28 $\xrightarrow{4.6 \text{ min}}$ 42:58

$\xrightarrow{5.0 \text{ min}}$ 42:58 $\xrightarrow{6.1 \text{ min}}$ 23:77 $\xrightarrow{6.15 \text{ min}}$ 0:100 $\xrightarrow{7.0 \text{ min}}$ 0:100 $\xrightarrow{7.05 \text{ min}}$ 90:10

Flow rate: 0.6 mL/min

Column: X-PressPak V-C18 (2- μm dia., 2.0 mm I.D. \times 50 mm)

Column temperature: 40
 Detection: Fluorescence; Ex 345 nm, Em 455 nm; Gain × 100


Figure 5. Analysis of standard mixture of 17 amino acids (20 pmol/μL each)
 This analysis was performed under the optimized conditions. Broken line indicates gradient profile (B%).


Figure 6. Linear dynamic ranges for 7 amino acid derivatives
 This figure indicates linear dynamic range for arbitrary chosen amino acid derivatives. Correlation coefficients (r) were also calculated for 7 amino acid derivatives ranging from 0.2 to 20 pmol. Excellent linearity was obtained.

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Table 1. Detection limit for 18 amino acid derivatives

| Amino acid | S/N=3 | Amino acid | S/N=3 |
|---------------------|---------|-----------------|---------|
| | fmol/μL | | fmol/μL |
| Cys-SO ₃ | 48.7 | Tyr | 20.9 |
| Asp | 40.9 | NH ₃ | 434.4 |
| Glu | 37.3 | Met | 23.6 |
| Ser | 30.5 | Val | 25.7 |
| His | 37.1 | Trp | 27.1 |
| Arg | 27.2 | Phe | 24.4 |
| Gly | 27.3 | Ile | 23.5 |
| Thr | 30.4 | Leu | 23.4 |
| Ala | 23.7 | Lys | 21.6 |

Highly sensitive detection was realized by the use of the fluorescence detector; the detection limit (S/N=3) was between 21 and 49 fmol/μL.


Figure 7. Analysis of amino acids in red wine
 This wine is a cooking wine.

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
Figure 8. Analysis of amino acids in functional beverage
This beverage contained 9 essential amino acids, i.e., His, Thr, Met, Val, Trp, Phe, Ile Leu and Lys.


Figure 9. Analysis of standard mixture of 20 amino acids containing GABA (20 pmol/μL each)
GABA (γ -aminobutyric acid) is an important inhibitory neurotransmitter in the central nervous system and essential for brain metabolism and function. Asparagine and glutamine are contained in this standard sample besides GABA. The analysis conditions are the same as shown in Fig. 5.

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