Technical Review

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Development of the Ultra Fast Liquid Chromatograph and the analytical column for high speed and high resolution analysis

Okiyuki Kunihiro, Yoshihiro Hayakawa, and Shuzo Maruyama

LC Business Unit, Analytical & Measuring Instruments Division, Shimadzu Corporation

1, Nishinokyo-Kuwabaracho Nakagyo-ku, Kyoto, 604-8511, Japan

Abstract

We developed Prominence UFLC as Ultra Fast Liquid Chromatograph and Shim-pack XR-ODS as ultra fast analytical column necessary to achieve higher speed and higher separation by using optimized particle size 2.2 μm . Prominence UFLC with Shim-pack XR-ODS columns enables users to shorten analysis time drastically and easily without extremely high pressure, while maintaining the high separation efficiency of conventional columns and system performance features such as reproducibility, carryover and durability. Prominence UFLC is applicable to a wide range of customers as a highly reliable, and flexible, fast LC system that does not require specialized components.

Keywords: ultra fast LC, optimized particle size, high speed and high separation, sample diffusion outside analytical column

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Fax : 075-811-8187

E-mail: kunihiro@shimadzu.co.jp


Figure 1. Relationship of linear velocity of mobile phase and HETP at each particle size


Figure 2. Pressure loss of the columns at each particle size

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Table 1. Theoretical plate number achieved with 30 MPa conditions

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
Figure 3. Relationship of retention time and particle size necessary to achieve $R_s=2$ under 30 MPa

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
Figure 4. Chromatogram of alkylphenones using Shim-pack XR-ODS (upper) and Shim-pack VP-ODS (lower) (Chromatographic Conditions : mobile phase: Water/Acetonitrile (3/7, v/v), temperature: 40 °C, detection: 245 nm, peak#: 1. acetophenone, 2. propiophenone, 3. butyrophenone, 4. valerophenone, 5. hexanophenone, 6. heptanophenone, 7. octanophenone)

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Table 2. Reproducibility of retention time and peak area at each peak (n=6)

Compounds	Retention time		Peak Area	
	Average	%RSD	Average	%RSD
Acetophenone	0.377	0.076	49018	0.193
Propiophenone	0.493	0.072	47260	0.167
Butyrophenone	0.587	0.070	49199	0.055
Valerophenone	0.677	0.065	45277	0.138
Hexanophenone	0.763	0.082	46613	0.186
Heptanophenone	0.842	0.087	43975	0.205
Octanophenone	0.914	0.080	48978	0.200

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
Figure 6. Chromatogram of alkylphenones after 100,000 injections
(Chromatographic condition: mobile phase: water/acetonitrile (1/1, v/v), flow rate: 1.25 mL/min, temperature: 40 °C Detection: absorbance at 245 nm)

Table 3. Reproducibility of retention time and peak area at each peak after 100,000 injections (n=6)

Compounds	Retention time		Peak Area	
	Average	%RSD	Average	%RSD
1. Acetophenone	0.916	0.089	312670	0.059
2. Propiophenone	1.398	0.058	315739	0.045
3. Butyrophenone	2.111	0.030	336428	0.063
4. Valerophenone	3.341	0.031	296609	0.040


Figure 7. Chromatogram of PTC-amino acids derivatives using Shim-pack XR-ODS
(Chromatographic Conditions : mobile phase: 10 mmol/L (potassium) phosphate <pH 7.0> for A and acetonitrile for B (A/B: 95/5 for 0.3 min then 95/5 to 60/40 in 3.4 min), flow rate: 1.2 mL/min, temperature: 40 °C, detection: absorbance at 254 nm)

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