Technical Review

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Retention behavior of newly developed ODS columns suitable for measurement of hydrophilic compounds

Hiroyuki Yamasaki, Seiji Itoh and Yasutami Mitoma

Separation Center, Scientific Instrument Division, Tosoh Corporation

Abstract

Octadecylsilane-bonded silica (ODS) columns are the most widely used columns for research and quality control in the pharmaceutical and chemical industries. We examined the retention behavior of various compounds on 2 kinds of newly developed ODS columns. Both packings were exhaustively reacted with end-capping reagent. Both columns showed excellent asymmetry factor with Desipramine, basic compound. It was confirmed that the retention of neutral compounds on ODS columns increased with ODS ligand density. On the other hand, TSKgel ODS-100 V, a column with low C 18 ligand density followed by exhaustive end-capping, strongly retained polar compounds compared with TSKgel ODS-100 Z, a column with high C 18 ligand density followed by exhaustive end-capping. Longer retention for polar compounds was assumed to be the ability of such analytes to interact with water-rich layers near the silica surface.

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Table 1. Properties of stationary phase

Column	Particle size (μm)	Pore size (Angstrom)	Specific surface area (m^2/g)	Pore volume (mL/g)	Functionality	Carbon contents ^{*)} (%)	Bonding structure
TSKgel ODS 100 V	5	100	450	1.10	C 18	15	Monolayer
TSKgel ODS 100 Z	5	100	450	1.10	C 18	20	Monolayer


Fig.1 Relation between Log P and Log k'


Fig.2 Column lifetime at low pH

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

Fig.3 Influence of pH of mobile phase on k'

Conditions
 Column :TSKgel ODS-100V
 Eluent :50mmol/L P.B.(pH X)/MeOH (30/70)
 Flow rate;1.0mL/min
 Temp. :40°C


Fig.4 Influence of pH of mobile phase on As'f

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
Fig 5 Chromatogram of benzene and desipramine

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- 1) Lold R. Snyder, Joseph J. Krikland, Joseph L. Glajch,
PRACTICAL HPLC METHOD DEVELOPMENT
- 2) N. Nagae and T. Enami, Bunseki Kagaku, 49 (2000) 887
- 3) N. Nagae, T. Enami, and S. Doshi, LC/GC Volume 10 Oct.
2002