Technical Review

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LC / Linear Ion Trap–Time of Flight Mass Spectrometer for Protein Analysis Using Nano LC With Dual Exchange Gradient System

Izumi Ogata, Kimihiko Ishii, Eri Yamashita, Shinji Yoshioka and Yasushi Terui

Hitachi High–Technologies Corporation, 882, Ichige, Hitachinaka, Ibaraki, 312–8504, Japan

Abstract

“NanoFrontier LD”, a new LC/Linear Ion Trap (LIT)– Time of Flight (TOF) MS for protein analysis has developed by Hitachi High–Technologies Corp. One of its remarkable features is NanoLC with Dual Exchange Gradient System (DEGS). DEGS offers great reproducibility of nano–flow gradient elution. Connecting NanoLC and LIT–TOF MS which MSⁿ analysis with high mass accuracy is available, Information Based Acquisition (IBA) is also available. With IBA, the efficiency of protein identification by MS/MS improves.

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Dual Exchange Gradient System (DEGS)

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Information Based Acquisition (IBA)

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
Figure 2. Schematic diagram of Dual Exchange Gradient System

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
Figure 3. Total Ion Chromatogram of BSA tryptic digest at 50 nL/min.


Figure 4. Total Ion Chromatograms of E.coli digest at 200 nL/min.

Table 1. Reproducibility of peak retention times of BSA tryptic digest

	Retention time (minutes)			
	peak 1	peak 2	peak 3	peak 4
exp. 1	15.4	27.4	31.4	40.9
exp. 2	15.3	27.2	31.4	40.8
exp. 3	15.4	27.3	31.2	40.9
exp. 4	15.3	27.3	31.4	40.5
exp. 5	15.6	27.3	31.0	40.9
mean	15.4	27.3	31.3	40.8
std. dev	0.095	0.090	0.188	0.172
rsd (%)	0.618	0.331	0.601	0.421

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
Figure 5. MS/MS Chromatograms of E.coli digest at 200 nL/min.


Figure 6. Statistics showing identified proteins for data number

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